# SOUTHERN BRAZILIAN JOURNAL OF CHEMISTRY Vol. 2, No.2, 1994

# NIOBIUM AS A POTENTIOMETRIC SENSOR IN REDOX TITRATIONS WITH AND WITHOUT PASSIVATION BY AMMONIUM MOLYBDATE.

Claudete J. Valduga, Eunice Valduga, Martha Adaime and Nádia Viaro<sup>1</sup>

Chemistry Department, Federal University of Santa Maria Santa Maria - RS - Brazil - 97119-900

#### ABSTRACT

The Pt indicator electrode has been commonly used for redox potentiometric titrations. In this work the Nb indicator electrode is suggested as an alternative for the Pt indicator electrode for the redox potentiometric titrations of the system  $Fe(II)/K_2Cr_2O_7$ . The utilization of proposed electrode has been studied with and without passivation. The passivation is made by electrode immersion in oxidising anodic inhibitor of NH4M0O4. The result was an observed increase of the potentiometric leap from 30 mV to 60 mV, allowing a better visualization of the equivalence point. The same equivalence point was observed with the conventional and proposed electrode. With the Nb electrode an inverse potentiometric leap was also obtained.

### RESUMO

O eletrodo indicador de platina tem sido usado nas titulações potenciométricas de oxi-redução convencionalmente. Neste trabalho, o eletrodo indicador de nióbio é considerado como uma alternativa ao eletrodo indicador de platina nas titulações potenciométricas de oxiredução para o sistema  $Fe(II)/K_2Cr_2O_4$ . A utilização do eletrodo foi estudada com e sem passivação. A passivação foi feita com a imersão do eletrodo num inibidor oxidante anódico de  $NH_4MoO_4$ . Como consequência ocorreu um aumento do salto potenciométrico de 30 mV para 60 mV (ECS), permitindo uma melhor visualização do ponto de equivalência. Foi observado o mesmo ponto de equivalência para o eletrodo proposto em relação ao convencional. Além disso, com o eletrodo de nióbio foi obtido um salto potenciométrico invertido.

# Keywords: potentiometric sensor; redox titrations; niobium; passivation; ammonium molybdate

<sup>&</sup>lt;sup>1</sup> Present Address: UFSM (campus), caixa postal 5051, Santa Maria 97119-911, RS-Brazil.

Niobium as Potentiometric Sensor

#### INTRODUCTION

The indicator electrodes used in potentiometric titrations depend on the type of reaction that is being investigated. In an acid-base titration it may be a hydrogen electrode or any other electrode that responds to hydrogen ions. The equivalence point of the reaction is revealed by a sudden change of potential in the e.m.f. curve vs the volume of the titrant solution.

Many pure metals are resistant to corrosion and have the common characteristic of forming oxides very adherent to the surface that repress the corrosion and are suitable as proton sensors<sup>1</sup>.

The metal-oxides. $xH_2O$  indicator electrodes have been used as sensors because of their low price and easy handling.

Some metals, such as Ti, Ni, Ta and alloys such as chromium-cobalt and stainless steel tend to become passive due to the formation of a thin, adherent layer of oxide or another compound that proctects them from oxidation<sup>2</sup>.

The sensibility of metal-oxide.xH<sub>2</sub>O electrodes in redox titration was studied by Sayed<sup>3</sup>. Redox titrations were made and appreciable potential leaps were observed in the equivalence point.

The Nb is covered by a Nb<sub>2</sub>O<sub>5</sub> film that can be dissolved by fusion with hydroxide but is not attacked by the hydrogen ion<sup>4</sup>.

Chromium-Nickel alloys, used as potentiometric sensor in redox titrations, show well determined potentiometric leaps, and usually in the same potential range of platinum electrodes. These potentiometric leaps are usually 50% of the value of the leaps obtained when platinum electrode are used<sup>5</sup>.

The first studies carried out to verify the response of metallic-oxide films used hastelloy alloys and have shown that the film can respond to protons in solution. Systematic studies in acid-base reactions and the Nernst type behavior extends to the range of pH 1,0 to pH 13,0<sup>6</sup>.

In this work, Nb was used as indicator electrode in redox potentiometric titration for the system  $Fe(II)/K_2Cr_2O_7$ . The passivation was made to increase the potentiometric leap.

# **MATERIALS AND METHODS**

Potentiometric titrations of 0,01N FeSO<sub>4</sub>.7H<sub>2</sub>O (in H<sub>2</sub>SO<sub>4</sub> solution) with 0,1N K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution were done, with stirring at room temperature using the calomel electrode as reference and a Nb electrode as indicator, instead of the Pt electrode. The Nb electrode was sanded with granulometry 280, 320, 400 and 600 a paper and then polished with Al<sub>2</sub>O<sub>3</sub>.

In some cases, the Nb electrodes were treated chemically and the passivation was achieved by immersion in 0,01 N NH4MoO4 solution for thirty minutes and for one hour.

#### SOUTH. BRAZ. J. CHEM., Vol. 2, No.2, 1994.

C.J. Valduga, E. Valduga. M. Adaime & N. Viaro

The titrations were conducted at  $25\pm$  1°C, with a PW 9409 Phillips potentiometer using a Shott Mainz B281 calomel electrode of simple junction, as reference. In order to achieve a better resolution a 274 microburete Methrohm was used. A 702 Fisaton magnetic stirrer, was used to ensure adequate mixing.

### **RESULTS AND DISCUSSION**

Figures 2 and 3 show that the equivalence points in the titrations coincide exactly with that verified when the Pt electrode is used (Figure 1), although the potential leap is reversed and smaller.

The Nb electrode passivated by immersion in oxidizing anodic inhibitor solution of NH<sub>4</sub>MoO<sub>4</sub> for thirty minutes and one hour presented a potentiometric leap of about 60 mV, that is twice that observed without passivation in the Fe(II) titration by  $K_2Cr_2O_7$  (30 mV). The passivated electrode allows a better detection of the equivalence point.

When the Nb is exposed to the air, it has the characteristic of formation of a film mainly of  $Nb_2O_5$ .

The oxidizing anodic inhibitor of NH<sub>4</sub>MoO<sub>4</sub>, probably makes this layer thicker and the redox potential becomes higher than the original one.




Figure 1. Titration curve using Pt indicator electrode for the system Fe (II)/Cr (VI).

57

Niobium as Potentiometric Sensor




Figure 2. Titration curve using Nb indicator electrode without passivation for the system Fe(II)/Cr (VI).




Figure 3. Titration curve using Nb indicator electrode with passivation for the system Fe (II)/Cr (VI).

#### SOUTIL BRAZ. J. CHEM., Vol. 2, No.2, 1994.

C.J. Valduga, E. Valduga, M. Adaime & N. Viaro

# REFERENCES

- 1. M. Ashrof-Chorassani and R.D. Braum, Corrosion NACE, 43, 32(1987).
- 2. J. Breme, "Titanium and Titanium Alloys the Biomaterial of Preference" (Plenary Lecture), 6th International Conference on Titanium, Cannes (1988).
- 3. S.M. Sayed, J. Chem., 21, 431 (1978).
- 4. R. C. Leite, "Nióbio, uma Conquista Nacional", 1ª Edição, Livraria Duas Cidades Ltda, São Paulo, S.P., 1988
- 5. A.C. Magalhães, M.D. Capelato, Uso de Ligas Ni-Cr como Sensor Potenciométrico em Titulações Redox, *Química Nova 14* (4), 37 (1991).
- 6. K. Nomura, Y.Ujhara, Anal. Chem., 60, 2564 (1988).

The SOUTHERN BRAZILIAN JOURNAL OF CHEMISTRY (ISSN: 2674-6891; 0104-5431) is an open-access journal since 1993. Journal DOI: 10.48141/SBJCHEM. http://www.sbjchem.com. This text was introduced in this file in 2021 for compliance reasons. © The Author(s)

OPEN ACCESS. This article is licensed under a Creative Commons Attribution 4.0 (CC BY 4.0) International License, which permits use, sharing, adaptation, distribution, and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third-party material in this article are included in the article's Creative Commons license unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit http://creativecommons.org/licenses/by/4.0/.